

## ORIHO CARE

# DAKSH MULTI TIBIA

**NAIL SYSTEM** 




**{ STANDARD QUALITY }** 





#### **CARE DAKSH Multi Tibia Nail Implants**



**DAKSH Multi Tibia Nail**Dia 8mm to 11mm Length 285mm to 390mm





#### **CARE DAKSH Multi Tibia Nail Instruments**



CDMT-1800.01 **Spanner** CDMT-1800.02 **Hexagonal Screw Driver Tip 3.5mm** CDMT-1800.03 Drill Bit ø 3.2 mm CDMT-1800.04 Drill Bit ø 4.0 mm **CDMT-1800.05 Protection Sleeve 10.0/8.0** CDMT-1800.06 **L-Key For Conical Bolt** CDMT-1800.07 Drill Sleeve 8.0/4.0 CDMT-1800.08 Trocar ø 8.0 mm CDMT-1800.09 **Aiming Device for Daksh Tibial Nail Along With Insertion Handle** 

Specifications, size, shape, color subject to change without notice. Picture shape if similar to any other manufacturers is just a co-incident.

CDMT-1800.10

**Guide RAM Cannulated** 





### **CARE DAKSH Multi Tibia Nail Instruments**



CDMT-1800.11 Depth Guage

CDMT-1800.12 Proximal Reamer

CDMT-1800.13 Driving Cap with thread, for

**Insertion Handle** 

CDMT-1800.14 Awl

CDMT-1800.15 Guide wire plain

CDMT-1800.16 Guide wire beaded

CDMT-1800.17 Tissue Protector

CDMT-1800.18 Medullary Tube





#### **CARE DAKSH Multi Tibia Nail Instruments**

CDMT-1800.19 Wire Holding Forcep

CDMT-18<mark>00.20 Hexagonal Screw Driver for</mark>

**Conical Bolt** 

CDMT-1800.21 Conical Bolt

Alliumine Alliumine

CDMT-1800.22 Measuring Scale

CDMT-1800.23 Slotted Hammer

CDMT-1800.24 Hammer



CDMT-1800.25 Sleeve for Proximal Reamer



**Container** 








CDMT - 1800.00

**DAKSH Multi Tibia Nail Instrument Set (Standard Quality)** 











# <u>CARE DAKSH Multi Tibia Nail Instrument Set</u> <u>(Standard Quality)</u>

Ref. No.	Description	Material	1 SET
CDMT-1800.00	DAKSH Multi Tibia Instrument Set (Standard Quality)		
	CONTENTS:-		
CDMT-1800.01	Spanner	SS	1
CDMT-1800.02	Hexagonal Screw Driver Tip 3.5mm	SS	1
CDMT-1800.03	Drill Bit ø 3.2 mm	SS	1
CDMT-1800.04	Drill Bit ø 4.0 mm	SS	1
CDMT-1800.05	Protection Sleeve 10.0/8.0,	SS	1
CDMT-1800.06	L-Key For Conical Bolt	SS	1
CDMT-1800.07	Drill Sleeve 8.0/4.0	SS	1
CDMT-1800.08	Trocar ø 8.0 mm	SS	1
CDMT-1800.09	Aiming Device for CARE Daksh Tibial Nail Along With	SS	1
	Insertion Handle		
CDMT-1800.10	Guide Ram Cannulated	SS	1
CDMT-1800.11	Depth Gauge	SS	1
CDMT-1800.12	Proximal Reamer	SS	1
CDMT-1800.13	Driving Cap with thread, for Insertion Handle	SS	1
CDMT-1800.14	Awl	SS	1
CDMT-1800.15	Guide wire plain	SS	1
CDMT-1800.16	Guide wire beaded	SS	1
CDMT-1800.17	Tissue Protector	SS	1
CDMT-1800.18	Medullary Tube	SS	1
CDMT-1800.19	Wire Holding Forcep	SS	1
CDMT-1800.20	Hexagonal Screw Driver for Conical Bolt	SS	1
CDMT-1800.21	Conical Bolt	SS	1
CDMT-1800.22	Measuring Scale	SS	1
CDMT-1800.23	Slotted Hammer	SS	1
CDMT-1800.24	Hammer	SS	1
CDMT-1800.25	Sleeve for Proximal Reamer	SS	1
	Container		

### www.orthocare.in www.orthoimplants.com





## ORIHO CARE

H.O.: 3/10 Kirti Nagar Industrial Area, New Delhi - 110015, INDIA Phone: +91-9810038246 | +91-11-41424511 | 41424522 E-mail: orthocareindia2009@gmail.com | Web: www.orthocare.in Factory: A-31, Sector A-4, Tronica City Loni, Ghaziabad (U.P.) INDIA

**USFDA 510K CLEARED PRODUCTS** CE 1023 CERIFIED PRODUCTS ISO 13485:2016 CERTIFICATION LICENSE NO. MFG/MD/2021/000347