



CAT NO. 659
 SINGLE PIN CLAMP



CAT NO. 658
 DOUBLE PIN CLAMP



CAT NO. 661
 OPEN CLAMP



CAT NO. 661M
 OPEN ADJUSTABLE
 CLAMP MODIFIED



CAT NO. 662
 CONNECTING CLAMP



CAT NO. 662A.100 TO 662A.300
 CONNECTING ROD 5mm
 Length : 100-150-200-250-300mm



CAT NO. 663
 TWIN ADJUSTABLE
 CLAMP-STRAIGHT



CAT NO. 663M
 OPEN ADJUSTABLE
 CLAMP MODIFIED STRAIGHT



CAT NO. 664
 TWIN ADJUSTABLE
 CLAMP-CURVED



CAT NO. 664 M
 OPEN ADJUSTABLE
 CLAMP MODIFIED CURVED



CAT NO. 664A
 TRANSVERSE CLAMP



CAT NO. 665
 TUBE TO TUBE CLAMP



CAT NO. 666
 SPANNER 11mm



CAT NO. 667
 T-WRENCH 11mm



CAT NO. 668
 COMPRESSION/DISTRACTION
 DEVICE



Best For The Life Time



1023

An ISO
 13485:2003
 WHO GMP



Tubular External Fixator



CAT NO. 669
SCHANZ SCREW
INTRODUCER



CAT NO. 670
TRIPLE TROCAR
WITH HANDLE



**CAT. NO. : CBF-515.11100 TO
CBF-515.11300**
CARBON : ROD
DIA : 11mm
LENGTH : 100mm TO 300mm



**CAT. NO. : CBF-515.11350 TO
CBF-515.11500**
CARBON : ROD
DIA : 11mm
LENGTH : 350mm TO 500mm



CAT NO. 671.100 TO 671.300
TUBULAR ROD Dia 11mm
Length : 100, 300mm



CAT NO. 671.350 to 671.600
TUBULAR ROD Dia 11mm
Length : 350, 500mm



CAT NO. 672
CIRCULAR TUBULAR ROD
Dia 11mm

Schanz Screw - Trocar Points



CAT NO. 676
SPANNER 8mm

CAT. NO.	DIA.	THREAD	LENGTH	CAT. NO.	DIA.	THREAD	LENGTH
673.01	4.5mm	25mm	150mm	673.21	4.5mm	50mm	250mm
673.02	4.5mm	25mm	175mm	673.22	5.0mm	50mm	125mm
673.03	4.5mm	25mm	200mm	673.23	5.0mm	25mm	150mm
673.04	4.5mm	50mm	225mm	673.24	5.0mm	50mm	150mm
673.05	5.0mm	25mm	200mm	673.25	5.0mm	25mm	175mm
673.06	5.0mm	50mm	225mm	673.26	5.0mm	50mm	175mm
673.07	6.0mm	50mm	200mm	673.27	5.0mm	40mm	180mm
673.08	4.0mm	25mm	150mm	673.28	5.0mm	32mm	200mm
673.09	4.0mm	25mm	175mm	673.29	5.0mm	50mm	200mm
673.10	4.0mm	32mm	175mm	673.30	5.0mm	25mm	225mm
673.11	4.0mm	25mm	200mm	673.31	5.0mm	50mm	225mm
673.12	4.5mm	32mm	175mm	673.32	5.5mm	25mm	175mm
673.13	4.5mm	40mm	175mm	673.33	5.5mm	25mm	200mm
673.14	4.5mm	50mm	175mm	673.34	5.5mm	50mm	200mm
673.15	4.5mm	32mm	200mm	673.35	5.0mm	60mm	200mm
673.16	4.5mm	50mm	200mm	673.36	5.0mm	25mm	170mm
673.17	4.5mm	25mm	225mm	673.37	5.0mm	40mm	175mm
673.18	4.5mm	32mm	225mm	673.38	5.5mm	60mm	200mm
673.19	4.5mm	32mm	250mm	673.39	5.0mm	25mm	250mm
673.20	4.5mm	50mm	225mm	673.40	6.0mm	50mm	100mm



Care For The Life Time



An ISO
13485:2003
WHO GMP





CAT NO. 674
CLAMPS FOR METCRAPLES
4mmX2.5mm



CAT NO. 678
SMALL COMPRESSION/
DISTRACTION DEVICE



CAT NO. 675
CLAMPS FOR FOREARM
4mmX3.5mm



CAT NO. 679
CONNECTING ROD 4mm
679.01 TO 679.05
Length: 100, 150, 200, 250, 300mm



CAT NO. 677
T-WRENCH 8mm

Schanz Screw - Trocar Points

CAT. NO.	DIA.	THREAD	LENGTH
680.00	2.0mm	12mm	100mm
680.01	2.5mm	25mm	100mm
680.02	3.5mm	25mm	100mm
680.03	3.0mm	18mm	100mm
680.04	3.0mm	25mm	125mm
680.05	3.5mm	18mm	100mm
680.06	3.5mm	25mm	100mm
680.07	3.5mm	32mm	125mm
680.08	3.5mm	32mm	150mm
680.09	1.5mm	12mm	100mm
680.10	3.0mm	18mm	75mm
680.11	3.0mm	25mm	110mm



Best For The Life Time



1023

An ISO
13485:2003
WHO GMP

